

Silane Compounds		
Product Code	Product Name	CAS No.
0101188-100G	AZIDOTRIMETHYL SILANE >97% puriss for synthesis (Trimethylsilylazide)	4648-54-8
0101188-25G	AZIDOTRIMETHYL SILANE >97% puriss for synthesis (Trimethylsilylazide)	4648-54-8
0101188-500G	AZIDOTRIMETHYL SILANE >97% puriss for synthesis (Trimethylsilylazide)	4648-54-8
0101193-100G	3-AMINOPROPYLTRIEHOXY SILANE (3-(Triethoxysilyl)-propylamine) >98% puriss for synthesis	919-30-2
0101193-500G	3-AMINOPROPYLTRIEHOXY SILANE (3-(Triethoxysilyl)-propylamine) >98% puriss for synthesis	919-30-2
0101317-10G	ALLYLTRIMETHYLSILANE >98% puriss for synthesis [762-72-1]	762-72-1
0101317-50G	ALLYLTRIMETHYLSILANE >98% puriss for synthesis [762-72-1]	762-72-1
0102127-100ML	N,O-BIS(TRIMETHYLSILYL) TRIFLUOROACETAMIDE(BSTFA) >98% puriss AR (Powerful silylating agent)	25561-30-2
0102127-25ML	N,O-BIS(TRIMETHYLSILYL) TRIFLUOROACETAMIDE(BSTFA) >98% puriss AR (Powerful silylating agent)	25561-30-2
0102127-5ML	N,O-BIS(TRIMETHYLSILYL) TRIFLUOROACETAMIDE(BSTFA) >98% puriss AR (Powerful silylating agent)	25561-30-2
0102169-100G	tert-BUTYLDIMETHYL SILYLCHLORIDE	18162-48-6
0102169-25G	tert-BUTYLDIMETHYL SILYLCHLORIDE	18162-48-6
0102169-500G	tert-BUTYLDIMETHYL SILYLCHLORIDE	18162-48-6
0102208-1000G	TERT-BUTYLCHLORODIPHENYLSILANE (tert-Butyldiphenylsilyl chloride) > 98% puriss for synthesis	58479-61-1
0102208-10G	TERT-BUTYLCHLORODIPHENYLSILANE (tert-Butyldiphenylsilyl chloride) > 98% puriss for synthesis	58479-61-1
0102208-250G	TERT-BUTYLCHLORODIPHENYLSILANE (tert-Butyldiphenylsilyl chloride) > 98% puriss for synthesis	58479-61-1
0102208-50G	TERT-BUTYLCHLORODIPHENYLSILANE (tert-Butyldiphenylsilyl chloride) > 98% puriss for synthesis	58479-61-1
0102218-100G	N,O-BIS(TRIMETHYLSILYL) ACETAMIDE (BSA) ~95% puriss for synthesis	10416-59-8
0102218-25G	N,O-BIS(TRIMETHYLSILYL) ACETAMIDE (BSA) ~95% puriss for synthesis	10416-59-8
0102218-500G	N,O-BIS(TRIMETHYLSILYL) ACETAMIDE (BSA) ~95% puriss for synthesis	10416-59-8
0102259-100G	BROMOTRIMETHYLSILANE > 98% puriss for synthesis [2857-97-8]	2857-97-8
0102259-25G	BROMOTRIMETHYLSILANE > 98% puriss for synthesis [2857-97-8]	2857-97-8
0102405-250G	BIS(TRIMETHYLSILYL)UREA >95% puriss for synthesis [18297-63-7]	18297-63-7
0102405-50G	BIS(TRIMETHYLSILYL)UREA >95% puriss for synthesis [18297-63-7]	18297-63-7
0102406-10G	BIS(TRIMETHYLSILYL)ACETYLENE >99% puriss for synthesis [14630-40-1]	14630-40-1
0102406-50G	BIS(TRIMETHYLSILYL)ACETYLENE >99% puriss for synthesis [14630-40-1]	14630-40-1
0102411-100G	TERT-BUTYLDIMETHYLSILANE >95% puriss for synthesis [29681-57-0]	29681-57-0
0102411-25G	TERT-BUTYLDIMETHYLSILANE >95% puriss for synthesis [29681-57-0]	29681-57-0
0102411-5G	TERT-BUTYLDIMETHYLSILANE >95% puriss for synthesis [29681-57-0]	29681-57-0
0102511-1G	1-(TERT-BUTYLDIMETHYLSILYL)IMIDAZOLE puriss for synthesis [54925-64-3]	54925-64-3
0102511-5G	1-(TERT-BUTYLDIMETHYLSILYL)IMIDAZOLE puriss for synthesis [54925-64-3]	54925-64-3
0102556-1G	(TERT-BUTYLDIMETHYLSILYL)ACETYLENE >98% puriss for synthesis [86318-61-8]	86318-61-8
0102556-25G	(TERT-BUTYLDIMETHYLSILYL)ACETYLENE >98% puriss for synthesis [86318-61-8]	86318-61-8
0102556-5G	(TERT-BUTYLDIMETHYLSILYL)ACETYLENE >98% puriss for synthesis [86318-61-8]	86318-61-8
0103306-100G	Chloromethyl trimethyl silane > 98% puriss for synthesis [2344-80-1]	2344-80-1
0103306-25G	Chloromethyl trimethyl silane > 98% puriss for synthesis [2344-80-1]	2344-80-1
0103306-500G	Chloromethyl trimethyl silane > 98% puriss for synthesis [2344-80-1]	2344-80-1
0104355-100ML	DIETHOXYDIMETHYLSILANE > 97% puriss for synthesis [78-62-6]	78-62-6
0104355-500ML	DIETHOXYDIMETHYLSILANE > 97% puriss for synthesis [78-62-6]	78-62-6
010441-1000ML	DIMETHYLDICHLOROSILANE 98% - 99%	75-78-5
010441-250ML	DIMETHYLDICHLOROSILANE 98% - 99%	75-78-5
0104420-100G	DIMETHYLPHENYLSILYLCHLORIDE >98% puriss for synthesis [768-33-2]	768-33-2
0104420-25G	DIMETHYLPHENYLSILYLCHLORIDE >98% puriss for synthesis [768-33-2]	768-33-2
0104420-5G	DIMETHYLPHENYLSILYLCHLORIDE >98% puriss for synthesis [768-33-2]	768-33-2
0104465-100G	N,N-DIETHYLAMINOTRIMETHYLSILANE >97% N-(Trimethylsilyl)diethylamine puriss for synthesis	996-50-9
0104465-25G	N,N-DIETHYLAMINOTRIMETHYLSILANE >97% N-(Trimethylsilyl)diethylamine puriss for synthesis	996-50-9
0104504-100G	DIMETHYLTHEXYL SILYL CHLORIDE puriss for synthesis [67373-56-2]	67373-56-2
0104504-25G	DIMETHYLTHEXYL SILYL CHLORIDE puriss for synthesis [67373-56-2]	67373-56-2
0105112-100ML	Ethyltrimethoxysilane > 97% puriss for synthesis [5314-55-6]	5314-55-6
0105112-500ML	Ethyltrimethoxysilane > 97% puriss for synthesis [5314-55-6]	5314-55-6
0106144-500G	FLUOROSILICIC ACID 30% puriss for synthesis [16961-83-4]	16961-83-4
010765-1000ML	3-(GLYCIDYLOXYPROPYL)TRIMETHOXSILANE >98% puriss for synthesis [2530-83-8]	2530-83-8
010765-250ML	3-(GLYCIDYLOXYPROPYL)TRIMETHOXSILANE >98% puriss for synthesis [2530-83-8]	2530-83-8

010810-100ML	1,1,1,3,3,3-HEXAMETHYLDISILAZANE (HMDS) >98% (GC)	999-97-3
010810-2500ML	1,1,1,3,3,3-HEXAMETHYLDISILAZANE (HMDS) >98% (GC)	999-97-3
010810-500ML	1,1,1,3,3,3-HEXAMETHYLDISILAZANE (HMDS) >98% (GC)	999-97-3
0108120-250G	HEXAMETHYLDISILANE > 99% puriss for synthesis [1450-14-2]	1450-14-2
0108120-50G	HEXAMETHYLDISILANE > 99% puriss for synthesis [1450-14-2]	1450-14-2
0108162-100ML	HEXAMETHYLDISILOXANE >98% puriss for synthesis [107-46-0]	107-46-0
0108162-500ML	HEXAMETHYLDISILOXANE >98% puriss for synthesis [107-46-0]	107-46-0
011233-100ML	LITHIUM-HMDS 20% (~1M) in THF (Lithium(trimethylsilyl)-amide), (Lithiumhexamethyldisilazide puriss	4039-32-1
011233-10x100ML	LITHIUM-HMDS 20% (~1M) in THF (Lithium(trimethylsilyl)-amide), (Lithiumhexamethyldisilazide puriss	4039-32-1
0113127-1000ML	METHYL TRICHLOROSILANE	75-79-6
0113127-250ML	METHYL TRICHLOROSILANE	75-79-6
0113181-100G	N-METHYL-N-(TRIMETHYLSILYL) TRIFLUOROACETAMIDE (MSTFA) >98% Derivatising agent used in GC analysis.	24589-78-4
0113181-25G	N-METHYL-N-(TRIMETHYLSILYL) TRIFLUOROACETAMIDE (MSTFA) >98% Derivatising agent used in GC analysis.	24589-78-4
0113181-5G	N-METHYL-N-(TRIMETHYLSILYL) TRIFLUOROACETAMIDE (MSTFA) >98% Derivatising agent used in GC analysis.	24589-78-4
0113216-500G	METHYLTRIS[METHYLETHYLKETOXIME]SILANE [MOS] [22984-54-9]	22984-54-9
0113250-100ML	Methyltriethoxysilane > 98% puriss for synthesis [2031-67-6]	2031-67-6
0113250-250ML	Methyltriethoxysilane > 98% puriss for synthesis [2031-67-6]	2031-67-6
0113251-1000ML	METHYLTRIMETHOXYSILANE > 97% puriss for synthesis [1185-55-3]	1185-55-3
0113251-250ML	METHYLTRIMETHOXYSILANE > 97% puriss for synthesis [1185-55-3]	1185-55-3
0113253-25G	N-(Methoxymethyl)-N-trimethylsilylmethyl benzylamine > 96% puriss for synthesis [93102-05-7]	93102-05-7
0113253-5G	N-(Methoxymethyl)-N-trimethylsilylmethyl benzylamine > 96% puriss for synthesis [93102-05-7]	93102-05-7
0113281-100ML	N-METHYL-N-(TERTBUTYLDIMETHYLSILYL)TRIFLUOROACETAMIDE (MBDSTFA) >99% [77377-52-7]	77377-52-7
0113281-25ML	N-METHYL-N-(TERTBUTYLDIMETHYLSILYL)TRIFLUOROACETAMIDE (MBDSTFA) >99% [77377-52-7]	77377-52-7
0113281-5ML	N-METHYL-N-(TERTBUTYLDIMETHYLSILYL)TRIFLUOROACETAMIDE (MBDSTFA) >99% [77377-52-7]	77377-52-7
0116257-1000G	Phenyltriethoxysilane > 98% puriss for synthesis [780-69-8]	780-69-8
0116257-250G	Phenyltriethoxysilane > 98% puriss for synthesis [780-69-8]	780-69-8
0116314-100ML	POTASSIUM HMDS 20 WT% IN THF (Potassium bis(trimethylsilyl)amide) puriss for synthesis [40949-94-8]	40949-94-8
0116314-500ML	POTASSIUM HMDS 20 WT% IN THF (Potassium bis(trimethylsilyl)amide) puriss for synthesis [40949-94-8]	40949-94-8
0116389-100G	POTASSIUM TRIMETHYL SILANOATE pract. for synthesis [10519-96-7]	10519-96-7
0116389-25G	POTASSIUM TRIMETHYL SILANOATE pract. for synthesis [10519-96-7]	10519-96-7
0116389-5G	POTASSIUM TRIMETHYL SILANOATE pract. for synthesis [10519-96-7]	10519-96-7
0116419-1G	PHENYL SILANE >98% puriss for synthesis [694-53-1]	694-53-1
0116419-25G	PHENYL SILANE >98% puriss for synthesis [694-53-1]	694-53-1
0116419-5G	PHENYL SILANE >98% puriss for synthesis [694-53-1]	694-53-1
0119211-10ML	SPECTROSIL 323 SILYLATING AGENT (equivalent of Sylon BTZ) for GC analysis	no CAS
011965-1000G	SPECTROSIL A (Silane coated colloidal silica) Density : 1.12	no CAS
011965-100G	SPECTROSIL A (Silane coated colloidal silica) Density : 1.12	no CAS
0120221-100G	TRIMETHYLSILYL IODIDE >97%	16029-98-4
0120221-25G	TRIMETHYLSILYL IODIDE >97%	16029-98-4
0120221-500G	TRIMETHYLSILYL IODIDE >97%	16029-98-4
0120261-100G	TRIETHYLCHLOROSILANE (TES-Chloride, Triethylsilylchloride) > 99% puriss for synthesis.	994-30-9
0120261-25G	TRIETHYLCHLOROSILANE (TES-Chloride, Triethylsilylchloride) > 99% puriss for synthesis.	994-30-9
0120261-5G	TRIETHYLCHLOROSILANE (TES-Chloride, Triethylsilylchloride) > 99% puriss for synthesis.	994-30-9
0120265-100G	TRIETHYLSILANE > 98% puriss for synthesis	617-86-7
0120265-25G	TRIETHYLSILANE > 98% puriss for synthesis	617-86-7
0120265-500G	TRIETHYLSILANE > 98% puriss for synthesis	617-86-7
0120313-100G	TRIISOPROPYLSILANE > 98% puriss for synthesis [6485-79-6]	6485-79-6
0120313-25G	TRIISOPROPYLSILANE > 98% puriss for synthesis [6485-79-6]	6485-79-6
0120313-5G	TRIISOPROPYLSILANE > 98% puriss for synthesis [6485-79-6]	6485-79-6
0120314-100G	TRIISOPROPYLSILYLCHLORIDE > 98% puriss for synthesis [13154-24-0]	13154-24-0
0120314-25G	TRIISOPROPYLSILYLCHLORIDE > 98% puriss for synthesis [13154-24-0]	13154-24-0
0120320-100G	TRIMETHYLSILYLACETYLENE > 98% puriss for synthesis [1066-54-2]	1066-54-2
0120320-25G	TRIMETHYLSILYLACETYLENE > 98% puriss for synthesis [1066-54-2]	1066-54-2
0120320-500G	TRIMETHYLSILYLACETYLENE > 98% puriss for synthesis [1066-54-2]	1066-54-2
0120333-100G	2-(Trimethylsilyl)ethanol > 99% puriss for synthesis [2916-68-9]	2916-68-9
0120333-25G	2-(Trimethylsilyl)ethanol > 99% puriss for synthesis [2916-68-9]	2916-68-9

0120333-5G	2-(Trimethylsilyl)ethanol > 99% puriss for synthesis [2916-68-9]	2916-68-9
0120334-1G	3,3,3-Trifluoropropyltrimethoxysilane > 97% puriss for synthesis [429-60-7]	429-60-7
0120334-5G	3,3,3-Trifluoropropyltrimethoxysilane > 97% puriss for synthesis [429-60-7]	429-60-7
0120338-100ML	N-Trimethylsilylimidazole > 98% puriss for synthesis [18156-74-6]	18156-74-6
0120338-25ML	N-Trimethylsilylimidazole > 98% puriss for synthesis [18156-74-6]	18156-74-6
0120338-5ML	N-Trimethylsilylimidazole > 98% puriss for synthesis [18156-74-6]	18156-74-6
0120356-100ML	Trimethylsilyldiazomethane 2M in Hexanes [18107-18-1]	18107-18-1
0120387-100G	(TRIFLUOROMETHYL)TRIMETHYLSILANE >99% puriss for synthesis [81290-20-2]	81290-20-2
0120387-25G	(TRIFLUOROMETHYL)TRIMETHYLSILANE >99% puriss for synthesis [81290-20-2]	81290-20-2
0120387-5G	(TRIFLUOROMETHYL)TRIMETHYLSILANE >99% puriss for synthesis [81290-20-2]	81290-20-2
0120392-25G	TRIBUTYLCHLOROSILANE (Chlorotributylsilane) >97% puriss for synthesis [995-45-9]	995-45-9
0120392-5G	TRIBUTYLCHLOROSILANE (Chlorotributylsilane) >97% puriss for synthesis [995-45-9]	995-45-9
0120393-10ML	TRIMETHYLSILANOL >98% puriss for synthesis [1066-40-6]	1066-40-6
0120393-50ML	TRIMETHYLSILANOL >98% puriss for synthesis [1066-40-6]	1066-40-6
0120394-100G	2-(TRIMETHYLSILYL)OXY-1,3-BUTADIENE puriss for synthesis [38053-91-7]	38053-91-7
0120394-25G	2-(TRIMETHYLSILYL)OXY-1,3-BUTADIENE puriss for synthesis [38053-91-7]	38053-91-7
0120394-5G	2-(TRIMETHYLSILYL)OXY-1,3-BUTADIENE puriss for synthesis [38053-91-7]	38053-91-7
0120427-1G	(TRIMETHYLSILYL)ACETYLENE >97% puriss for synthesis [1777-03-3]	1777-03-3
0120427-5G	(TRIMETHYLSILYL)ACETYLENE >97% puriss for synthesis [1777-03-3]	1777-03-3
0120428-1G	(TRIISOPROPYLSILYL)ACETYLENE >97% puriss for synthesis [89343-06-6]	89343-06-6
0120428-25G	(TRIISOPROPYLSILYL)ACETYLENE >97% puriss for synthesis [89343-06-6]	89343-06-6
0120428-5G	(TRIISOPROPYLSILYL)ACETYLENE >97% puriss for synthesis [89343-06-6]	89343-06-6
0120429-1G	1-(TRIMETHYLSILYL)-1-PROPYLENE >99% puriss for synthesis [6224-91-5]	6224-91-5
0120429-25G	1-(TRIMETHYLSILYL)-1-PROPYLENE >99% puriss for synthesis [6224-91-5]	6224-91-5
0120429-5G	1-(TRIMETHYLSILYL)-1-PROPYLENE >99% puriss for synthesis [6224-91-5]	6224-91-5
0120468-1G	1-TRIISOPROPYLSILYL-1-PROPYLENE puriss for synthesis [82192-57-2]	82192-57-2
0120468-5G	1-TRIISOPROPYLSILYL-1-PROPYLENE puriss for synthesis [82192-57-2]	82192-57-2
0120520-100G	1,1,3,3-TETRAMETHYLDISILOXANE >97% puriss for synthesis [3277-26-7]	3277-26-7
0120520-25G	1,1,3,3-TETRAMETHYLDISILOXANE >97% puriss for synthesis [3277-26-7]	3277-26-7
0120524-2500ML	TETRAETHYLORTHO SILICATE >99% puriss for synthesis [78-10-4]	78-10-4
0120524-500ML	TETRAETHYLORTHO SILICATE >99% puriss for synthesis [78-10-4]	78-10-4
0120528-100ML	1,1,3,3-TETRAMETHYLDISILOXANE >98% puriss for synthesis	3277-26-7
0120528-25ML	1,1,3,3-TETRAMETHYLDISILOXANE >98% puriss for synthesis	3277-26-7
012053-1000ML	TRIMETHYLCHLOROSILANE >98% (GC)	75-77-4
012053-250ML	TRIMETHYLCHLOROSILANE >98% (GC)	75-77-4
012057-100ML	TRIMETHYLSILYLCYANIDE >98%	7677-24-9
012057-25ML	TRIMETHYLSILYLCYANIDE >98%	7677-24-9
012212-100ML	Vinyltriethoxysilane > 98% puriss for synthesis [78-08-0]	78-08-0
012212-500ML	Vinyltriethoxysilane > 98% puriss for synthesis [78-08-0]	78-08-0